

Supporting Information

Rec. Nat. Prod. 14:5 (2020) 383-386

Scoparic acid E: A New Labdane Diterpenoid on Attenuating Palmitate Induced Viability in MIN6 cells from *Scoparia dulcis*

Caiyun Zhang^{1#}, Lizhen Chen^{1#}, Yiyi Li², Na Wei¹, Lu Zhang¹, Lin
Dong¹, Yong Zhang^{*2} and Xiaopo Zhang^{*1}

¹Key Laboratory of Tropical Translational Medicine of the Ministry of Education, Hainan Key
Laboratory for Research and Development of Tropical Herbs, School of Pharmacy, Hainan
Medical University, Haikou 571199, P. R. China

²Department of Basic Medicine and Life Sciences, Hainan Medicinal University, Haikou 571199,
P. R. China

Table of Contents	Page
Figure S1: HRESIMS Spectrum of 1 (Scoparic acid E)	2
Figure S2: ¹ H-NMR (600 MHz, CDCl ₃) Spectrum of 1 (Scoparic acid E)	3
Figure S3: Enlarged ¹ H-NMR (600 MHz, CDCl ₃) Spectrum of 1 (Scoparic acid E)	4
Figure S4: ¹³ C-NMR (150 MHz, CDCl ₃) Spectrum of 1 (Scoparic acid E)	4
Figure S5: Enlarged ¹³ C-NMR (150 MHz, CDCl ₃) Spectrum of 1 (Scoparic acid E)	5
Figure S6: Enlarged ¹³ C-NMR (150 MHz, CDCl ₃) Spectrum of 1 (Scoparic acid E)	5
Figure S7: HSQC Spectrum of 1 (Scoparic acid E)	6
Figure S8: Enlarged HSQC Spectrum of 1 (Scoparic acid E)	6
Figure S9: ¹ H- ¹ H COSY Spectrum of 1 (Scoparic acid E)	7
Figure S10: Enlarged ¹ H- ¹ H COSY Spectrum of 1 (Scoparic acid E)	7
Figure S11: HMBC Spectrum of 1 (Scoparic acid E)	8
Figure S12: Enlarged HMBC Spectrum of 1 (Scoparic acid E)	8
Figure S13: ROESY Spectrum of 1 (Scoparic acid E)	9
Figure S14: Scifinder Search Report for 1 (Scoparic acid E)	9
Figure S15: Scifinder Search Report of 95% Similar Compound of 1 (Scoparic acid E)	10
Table S1: NMR data of Scoparic acid E and C (δ in ppm, J values in Hz)	11
S1: Detail of Bioactivity Test-Cell Viability Assay	12

* Corresponding authors: e-mail: 935150242@qq.com (Y. Zhang); E-mail: z_xp1412@163.com (X. P. Zhang)

Authors contributed equally to this paper.




Figure S1: HRESIMS Spectrum of **1** (Scoparic acid E)




Figure S2: ¹H-NMR (600 MHz, CDCl₃) Spectrum of **1** (Scoparic acid E)




Figure S3: Enlarged $^1\text{H-NMR}$ (600 MHz, CDCl_3) Spectrum of **1** (Scoparic acid E)




Figure S4: $^{13}\text{C-NMR}$ (150 MHz, CDCl_3) spectrum of **1** (Scoparic acid E)




Figure S5: Enlarged ^{13}C -NMR (150 MHz, CDCl_3) spectrum of **1** (Scoparic acid E)




Figure S6: Enlarged ^{13}C -NMR (150 MHz, CDCl_3) spectrum of **1** (Scoparic acid E)




Figure S7: HSQC spectrum of 1 (Scoparic acid E)




Figure S8: Enlarged HSQC spectrum of 1 (Scoparic acid E)




Figure S9: ^1H - ^1H COSY spectrum of **1** (Scoparic acid E)




Figure S10: Enlarged ^1H - ^1H COSY spectrum of **1** (Scoparic acid E)




Figure S11: HMBC spectrum of **1** (Scoparic acid E)




Figure S12: Enlarged HMBC spectrum of **1** (Scoparic acid E)




Figure S13: ROESY Spectrum of **1** (Scoparic acid E)

Select All Deselect All

0 of 8 Similarity Candidates Selected

	Substances
<input type="checkbox"/> ≥ 99 (most similar)	0
<input type="checkbox"/> 95-98	1
<input type="checkbox"/> 90-94	4
<input type="checkbox"/> 85-89	8
<input type="checkbox"/> 80-84	6
<input type="checkbox"/> 75-79	162
<input type="checkbox"/> 70-74	1778
<input type="checkbox"/> 65-69	10823
<input type="checkbox"/> 0-64 (least similar)	52711

Get Substances

Figure S14: Scifinder Search Report for 1 (Scoparic acid E)




Figure S15: Scifinder Search Report for 95% Similar Compound 1 (Scoparic acid E)

Table S1: NMR Data of Scoparic acid E and Scoparic acid C (δ in ppm, J values in Hz)

No.	Scoparic acid E ^a		Scoparic acid C	
	δ_C	δ_H	δ_C	δ_H
1	37.8	1.74, m 1.20, m	38.16	
2	18.2	1.78, m 1.61, m	18.79	
3	40.0	1.72, m	40.21	
4	47.0		48.18	
5	42.3	2.65, brs	44.08	2.42, brs
6	74.0	5.47, brs	74.10	4.71, brs
7	37.0	2.60, brd, (13.2) 2.48, brd, (13.2)	37.61	1.97, m
8	144.4		144.40	
9	57.4	1.93, brs	57.71	
10	38.6		38.55	
11	26.8	2.10, m 2.22, t, (13.2)	24.28	
12	32.4	2.00, m 2.15, m	26.70	
13	140.7		150.39	
14	172.9		194.82	
15				
16	127.6	5.67, s 6.26, s	134.23	5.91, s 6.26, s
17	113.7	4.76, m	113.14	4.65, brs
18	184.7		185.30	
19	19.2	1.36, s	19.73	1.07, s
20	25.6	1.50, s	25.53	1.36, s
1'	129.6		130.73	
2', 6'	130.7	8.04, d, (8.4)	129.88	7.95, d, (7.3)
3', 5'	128.4	7.46, t, (8.4)	128.32	7.30, t, (7.3)
4'	132.8	7.56, t, (8.4)	132.89	7.46, t, (7.3)
7'	166.2		167.35	

^a Measured at 600 MHz for ¹H NMR and 150 MHz for ¹³C NMR in CDCl₃.

S1: Detail of Bioactivity Test-Cell Viability Assay

Cell Lines and Cell Culture

The MIN6 cells were cultured in DRPMI 1640 Medium equilibrated with 5% CO₂ and 95% air at 37°C. The medium was supplemented with 10% fetal calf serum, 100 U/ml penicillin sulfate and 50 µg/ml gentamycin. All experiments were performed when cells reached 80%-90% confluence.

Cell Viability Assay

Cell viability was assessed by the MTT assay. Briefly, MIN6 cells were seeded in 96 well plates at 1×10^4 cells/well. Cells were incubated with 300 µM palmitate for 24 h. MTT solution was added to the cells at a final concentration of 0.5 mg/ml. After incubating for 4 h at 37 °C, with 5% CO₂, the solution was removed, and 150 µL DMSO was added. The precipitate in each well was dissolved for 10min and the optical density (OD) was determined at 570 nm using a microplate reader. The viability of the new compound treated cells was performed as above. The cell viability was calculated according to the following formula.

$$\text{Cell viability (\%)} = \frac{\text{PA (OD)} - \text{Normal (OD)}}{\text{Compound (OD)} - \text{Normal (OD)}} \times 100\%$$