

Supporting Information

J. Chem. Metrol. 16:2 (2022) 111-124

Validated chromatographic methods for concurrent determination of atorvastatin and perindopril

Sagar Patel, Chandani Patel Jasmina Surati, Ashok Akabari
and Ketan Shah

Shree Naranjibhai Lalbhai Patel College of Pharmacy, Umrakh, Gujarat, India

Table of Content	Page No.
Figure S1: IR spectra of atorvastatin	2
Figure S2: IR spectra of perindopril	2
Figure S3: Photograph of calibration curve of ato (800-4800ng/band), per (200-1200 ng/band)	3
Figure S4: Blank chromatogram of mobile phase	3
Figure S5: Densitogram of atorvastatin (2400ng/band), perindopril (600ng/band) using dichloromethane: methanol: ethyl acetate: glacial acetic acid (6:2:2:0.1, v/v/v/v) as a mobile phase	4
Figure S6: Chromatogram of ATO (40µg/mL) and PER (10µg/mL) in acetonitrile: methanol: buffer (40:10:50v/v/v) (ph 3) having RT of ATO at 4.39 min and PER at 7.82 min	4
Figure S7: Peak purity spectra of calibration range	5
Figure S8: Chromatogram of ATO (10µg/mL) in acetonitrile: methanol: Buffer (40:10:50v/v/v) (pH 3)	5
Figure S9: Chromatogram of PER (10µg/mL) in acetonitrile: methanol: buffer (40:10:50v/v/v) (pH 3)	6
Figure S10: Chromatogram of mobile phase	6
Table S1: Linearity of atorvastatin and perindopril (HPTLC method)	7
Table S2: Result of intermediate precision of atorvastatin (HPTLC method)	7
Table S3: Result of intermediate precision of perindopril (HPTLC method)	7
Table S4: Result of repeatability precision of atorvastatin (HPTLC method)	8
Table S5: Result of repeatability precision of perindopril (HPTLC method)	8
Table S6: Accuracy data of atorvastatin (HPTLC method)	8
Table S7: Accuracy data of perindopril (HPTLC method)	8
Table S8: Assay result of synthetic mixture (HPTLC method)	9
Table S9: Linearity of atorvastatin and perindopril (HPLC method)	9
Table S10: Intraday precision data for estimation of ATO and PER (n=3) (HPLC method)	9
Table S11: Interday precision data for estimation of ATO and PER (n=3) (HPLC method)	9
Table S12: Repeatability data for ATO and PER (HPLC method)	10
Table S13: Accuracy data of ATO (HPLC method)	10
Table S14: Accuracy data of PER (HPLC method)	10
Table S15: Analysis data of formulation (HPLC method)	11
Table S16: F-Test two-sample for variances	11




Figure S1: IR Spectra of atorvastatin




Figure S2: IR Spectra of perindopril




Figure S3: Photograph of calibration curve of ATO A1-A6 (800-4800ng/band), PER P1-P6 (200-1200 ng/band)




Figure S4: Blank chromatogram of mobile phase




Figure S5: Densitogram of atorvastatin (2400ng/band), perindopril (600ng/band) using dichloromethane: methanol: ethyl acetate: glacial acetic acid (6:2:2:0.1, v/v/v/v) as a mobile phase




Figure S6: Chromatogram of ATO (40 μ g/mL) and PER (10 μ g/mL) in acetonitrile: methanol: buffer (40:10:50v/v/v) (ph 3) having RT of ATO at 4.39 min and PER at 7.82 min




Figure S7: Peak purity spectra of calibration range




Figure S8: Chromatogram of ATO (10µg/ml) in Acetonitrile: Methanol: Buffer(40:10:50v/v/v) (pH 3)




Figure S9: Chromatogram of PER (10µg/ml) in Acetonitrile: Methanol: Buffer(40:10:50v/v) (pH 3)




Figure S10: Chromatogram of mobile phase

Table S1: Linearity of atorvastatin and perindopril (HPTLC method)

Sr. no	Concentration (ng/band)	Perindopril		Atorvastatin	
		Mean peak area $\pm$ S.D. (n=6)	R.S.D. (%)	Mean peak area $\pm$ S.D. (n=6)	R.S.D. (%)
1	800	5899.2 $\pm$ 86.78	1.84	630.01 $\pm$ 9.98	1.73
2	1600	10432.7 $\pm$ 108.21	1.71	1368.21 $\pm$ 3.86	0.28
3	2400	15821.5 $\pm$ 89.62	0.57	2310.01 $\pm$ 21.90	0.95
4	3200	19978.7 $\pm$ 140.29	0.78	3038.63 $\pm$ 57.62	1.85
5	4000	24987.2 $\pm$ 155.84	1.13	4138.1 $\pm$ 51.42	1.23
6	4800	29958.2 $\pm$ 202.80	1.02	5002.35 $\pm$ 58.36	1.12

Table S2: Result of intermediate precision of atorvastatin (HPTLC method)

Atorvastatin Conc.(ng/band)	Intraday		Interday	
	Area $\pm$ SD (n=3)	%RSD	Area $\pm$ SD (n=3)	%RSD
800	5770.86 $\pm$ 91.62	1.58	5781.7 $\pm$ 92.42	1.21
2400	15847.83 $\pm$ 72.19	0.45	15645.43 $\pm$ 69.78	0.68
4800	30079.17 $\pm$ 423.27	1.40	30125.17 $\pm$ 356.04	1.24

Table S3: Result of intermediate precision of perindopril (HPTLC method)

Perindopril Conc. (ng/band)	Intraday		Interday	
	Area $\pm$ SD (n=3)	%RSD	Area $\pm$ SD (n=3)	%RSD
200	635.83 $\pm$ 4.64	0.73	629.13 $\pm$ 3.96	0.98
600	3073.03 $\pm$ 41.40	1.34	3056.41 $\pm$ 31.40	1.24
1200	5051.13 $\pm$ 31.77	0.62	5123.4 $\pm$ 29.76	1.65

Table S4: Result of repeatability precision of atorvastatin (HPTLC method)

Injection			Rf value	
Conc.ng/band	Area $\pm$ SD (n=6)	%RSD	Rf $\pm$ SD (n=6)	%RSD
2400	11055.65 $\pm$ 162.95	1.47	0.70 $\pm$ 0.0037	0.52

Table S5: Result of repeatability precision of perindopril (HPTLC method)

Injection			Rf value	
Conc.ng/band	Area $\pm$ SD(n=6)	%RSD	Rf $\pm$ SD (n=6)	%RSD
600	1414.55 $\pm$ 27.04	1.91	0.39 $\pm$ 0.0035	0.95

Table S6: Accuracy data of atorvastatin (HPTLC method)

Level	Conc. of ATO from Synthetic mixture (μ g/ml)	Amount of Std. ATO added (μ g/ml)	Total amount of ATO (μ g/ml)	Total amount of ATO Recovered (μ g/ml) * Mean $\pm$ SD	% Recovery	% RSD
50%	1600	800	2400	2428.4 $\pm$ 84.47	101.18%	0.028
100%	1600	1600	3200	32183.9 $\pm$ 415.3	100.59%	1.97
150%	1600	2400	4000	2607.9 $\pm$ 983.33	100.65%	1.58

Table S7: Accuracy data of perindopril (HPTLC method)

Level	Conc. of PER from Synthetic mixture (μ g/ml)	Amount of Std. PER added (μ g/ml)	Total amount of PER (μ g/ml)	Total amount of PER Recovered (μ g/ml) Mean $\pm$ SD	% Recovery	% RSD
50%	400	200	600	605.34 $\pm$ 9.52	99.79%	0.41
100%	400	400	800	801.27 $\pm$ 272.28	100.54%	1.21
150%	400	600	1000	1011 $\pm$ 7.37	101.1%	0.98

Table S8: Assay result of synthetic mixture (HPTLC method)

Formulation	Drug	Conc. (ng/band)	%Assay± SD (n=3)	% R.S.D
Synthetic	Atorvastatin	2400	101.6±30.15	1.21
mixture	Perindopril	600	100.42 ± 45.34	1.76

Table S9: Linearity of Atorvastatin and Perindopril (HPLC method)

Sr. no	Concentration (µg/ml)	Perindopril		Atorvastatin	
		Mean peak area ± S.D. (n=6)	R.S.D. (%)	Mean peak area ± S.D. (n=6)	R.S.D. (%)
1	10	5985±105.32	1.73	14986±289.98	1.90
2	20	12452±225.30	1.83	31246±328.88	1.05
3	30	17895±314.26	1.73	45095±351.54	0.76
4	40	24423±202.85	0.84	61562±390.03	0.63
5	50	29895±194.57	0.64	75426±400.45	0.52

Table S10: Intraday precision data for estimation of ATO and PER (n=3) (HPLC method)

Conc. (µg/ml)		Mean peak area	% RSD	Mean peak area	% RSD
		±SD		±SD	
ATO	PER	ATO		PER	
20	10	15063±269.22	1.78	6109.33±88.46	1.44
60	30	45076±98.04	0.21	18002.67±228.55	1.26
100	50	76025.33±136.10	0.17	28031.33±308.03	1.09

Table S11: Interday precision data for estimation of ATO and PER (n=3) (HPLC method)

Conc. (µg/ml)		Mean peak area	% RSD	Mean peak area ±	% RSD
		±SD		SD	
ATO	PER	ATO		PER	
20	10	15213±207.88	1.36	5929.33±14.07	0.23
60	30	45348.33±303.93	0.67	18017.33±221.69	1.23
100	50	176350±339.04	0.44	29246.67±236.51	0.80

Table S12: Repeatability data for ATO and PER (HPLC method)

Drug	Concentration (µg/ml)	Mean peak area ± S.D. (n=7)	% R.S.D.
ATO	40	31837.83±327.53	1.02
PER	20	12671.33±252.68	1.99

Table S13: Accuracy data of ATO (HPLC method)

Level	Conc. of ATO from Synthetic mixture (µg/ml)	Amount of Std. ATO added (µg/ml)	Total amount of ATO (µg/ml)	Total amount of ATO Recovered (µg/ml) * Mean ± SD	% Recovery	% RSD
50%	40	20	60	59.38±399.24	98.96%	0.87
100%	40	40	80	80.43±143.77	100.53%	0.23
150%	40	60	100	99.33±266.37	99.33%	0.35

Table S14: Accuracy data of PER (HPLC method)

Level	Conc. of PER from Synthetic Mixture (µg/ml)	Amount of Std. PER added (µg/ml)	Total amount of PER (µg/ml)	Total amount of PER Recovered (µg/ml) Mean ± SD	% Recovery	% RSD
50%	10	5	15	15.07±25.31	100.52%	0.27
100%	10	10	20	20.59±61.75	100.95%	0.49
150%	10	15	25	237.36±25.08	100.33%	1.56

Table S15: Analysis data of formulation (HPLC method)

Sr. No	Drug	Concentration (µg/ml)	% Assay ± SD	%R.S. D
1	ATO	40	100.98±200.31	0.60
2	PER	10	100.00 ±60.53	0.46

Table S16: F-test two-sample for variances

	Perindopril		Atorvastatin	
	HPLC	HPTLC	HPLC	HPTLC
Mean	100.446	100.15	100.12	100.3
Variance	0.18628	0.10885	0.467	0.425
Observations	5	5	5	5
df	4	4	4	4
F	1.711345889		1.098824	
P(F<=f) one-tail	0.307745467		0.464712	
F Critical one-tail	6.388232909		6.388233	