

Supporting Information

J. Chem. Metrol. X:X (202X) XX-XX

Application of the localized surface plasmon resonance of silver nanoparticles for the determination of ethion in water samples

Hesam Daneshyar ^{1,2}, Seyyed Hamid Ahmadi ^{2*}
and Naser Dalali ¹

¹*Phase Separation & FIA Lab., Department of Chemistry, Faculty of Science, University of Zanjan, Zanjan, Iran*

²*Department of Environmental Analytical Chemistry, Faculty of Clean Technologies, Chemistry and Chemical Engineering Research Center of Iran, Tehran, Iran*

Table of Contents	Page
Figure S1: Chemical structure of ethion.	2
Figure S2: (a) Typical absorption spectra of Ag-NPs (0.25 mM) in presence of different concentration of ET (1-40 μ M) (b) Calibration curve of LSPR signal for ET in optimum conditions.	3




Figure S1: Chemical structure of ethion.




Figure S2: (a) Typical absorption spectra of Ag-NPs (0.25 mM) in presence of different concentration of ET (1-40 μM) (b) Calibration curve of LSPR signal for ET in optimum conditions